

Background of the Invention

The present invention relates to the field of the disinfection and decontamination of equipment and
5 instruments having surfaces liable to act as support for the deposition of a biofilm.

Description of the Prior Art

10 This equipment and/or these instruments are, for example, those used in the medical field, such as analytical instruments and any equipment such as nonautoclavable reusable medical devices, such as
15 dialysis generators, and also implants (ocular implants, heart valves) and prostheses. The equipment used in dentistry, even mucus membranes and teeth themselves, whether natural or prosthetic, are also liable to be the site of biofilm deposition. The
20 equipment used in the food and/or pharmaceutical industry may also be cited, as well as air-conditioning plants and, more generally, any equipment in contact with an aqueous medium liable to contain bacteria in suspension.

25 The current major problem is the removal of the biofilm formed by the biomass attached to the surfaces of the equipment, instruments and/or mucus membranes, which biofilm constitutes a cause of persistent infections and/or of contamination. This is because any bacterium
30 in suspension in an aqueous medium has the property of adhering to the supports that it encounters in order to form a biofilm. This biofilm is an agglomerate of bacteria on a surface, these bacteria being surrounded by a matrix of exopolysaccharides; formation of the
35 biofilm is a natural phenomenon. Once the biofilm has been formed, it is very difficult to remove.

Current decontamination and/or disinfection methods that have been proposed for combating biofilms, although

having a certain degree of effectiveness especially from the antibacterial standpoint, do not, however, remove the biofilm from the support, thereby encouraging its redevelopment and, in the particular case of hemodialysis, leaves supports for pyrogens on the surface.

In certain instruments or equipment, the biofilm is also enhanced by deposits of scale formed from calcium carbonate or magnesium carbonate.

The current use of decalcifying solutions to enhance the action of purely disinfecting solutions for the treatment of certain equipment does not, however, completely remove this biofilm except by using products such as bleach, which products, although effective against the biofilm, can often destroy the medical instruments and equipment. Furthermore, these solutions cannot be used on surfaces in contact with tissue and/or directly on mucus membranes.

L.F. Jacquelin, Pathologie Biologie, May 1994, page 425, teaches the sequential use of enzymes and of a phenolic disinfectant for the destruction of biofilms.

C. Johansen, Applied and Environmental Microbiology, September 1997, page 3724, teaches the use of enzyme combinations such as glucose oxydases and lactoperoxydase. WO 01/53010 teaches an enzymatic method of removing biofilms, comprising the use of an enzyme belonging to the group of carbohydrases and proteases and their sequential use and also in combination with, or independently of, agents belonging to the group of biocides, chelating agents and other cleaning agents. Also taught, in EP 1 186 574, is a method of removing biofilms from surfaces in contact with water, characterized in that it includes cleaning with an enzymatic active principle and cleaning with a disinfectant intended to kill the bacteria released by the action of the enzyme mixture, but the results

obtained are unsatisfactory and, at the present time, there is no process or product able to remove these biofilms.

- 5 The present invention allows the problem to be solved by the implementation of a method and a selection of products making it possible to obtain a level of removal effectiveness never hitherto achieved.

10 Summary of the Invention

The present invention relates to a method of removing a biofilm, which comprises at least the following steps, carried out simultaneously or consecutively:

- 15 a) a solution comprising an enzyme mixture containing at least one enzyme chosen from the group of proteases and at least one enzyme chosen from the group of osidases is prepared;
- b) a solution comprising a detergent exhibiting
- 20 antibacterial, antiviral and antifungal properties is prepared; and
- c) said solutions are applied, by washing or by circulation, to the surface to be treated.

- 25 The method according to the invention is one which furthermore comprises the following steps, carried out simultaneously or consecutively:

- d) a solution comprising an acid capable of dissolving deposits of mineral salts is prepared; and
- 30 e) said solution is applied, by washing or by circulation, to the surface to be treated.

The invention also relates to:

- the method according to the invention, wherein
- 35 the protease is chosen from the group formed by exopeptidases or endopeptidases, such as trypsin;
- the method according to the invention, wherein the osidase or carbohydrase is chosen from the group formed by amylase, glycosidase and galactosidase;

- the method according to the invention, wherein the enzyme mixture furthermore comprises an enzyme chosen from the group formed by hydrolases that break ester bonds, or esterases, especially carboxyl ester
5 hydrolases, such as lipase, phospholipases and/or phosphonodiesterases, such as ribonuclease;

- the method according to the invention, wherein the enzyme mixture is pancreatin;

- the method according to the invention, wherein
10 the detergent exhibiting antibacterial, antiviral and antifungal properties is a solution containing surfactants;

- the method according to the invention, wherein the detergent exhibiting antibacterial, antiviral and
15 antifungal properties is an alkaline solution containing surfactants and a quaternary ammonium;

- the method according to the invention, wherein the detergent exhibiting antibacterial, antiviral and
20 antifungal properties is a neutral solution containing surfactants;

- the method according to the invention, wherein the detergent exhibiting antibacterial, antiviral and
antifungal properties is an acid solution containing surfactants, in this variant, since the detergent is an
25 acid solution, it possesses both antibacterial, antiviral and antifungal properties and allows the deposits of mineral salts to be dissolved without an additional step; and

- the method according to the invention, wherein
30 the detergent furthermore contains a disinfectant such as a sodium hypochlorite or potassium hypochlorite solution.

When, in a variant, the method comprises a step of
35 washing with an acid solution in order to remove the deposits of mineral salts, the acid is chosen from the group formed by citric acid, peractetic acid, glycolic acid and oxyacetic acid.

The invention also relates to a kit intended for removing a biofilm, which comprises at least one solution of an enzyme mixture containing at least one enzyme chosen from the group of proteases and at least one enzyme chosen from the group of osidases, and at least one solution of a detergent exhibiting antibacterial, antiviral and antifungal properties.

The invention also relates to:

- 10 - a kit according to the invention, wherein the protease is chosen from the group formed by exopeptidases or endopeptidases, such as trypsin;
- a kit according to the invention, wherein the osidase or carbohydrase is chosen from the group formed
- 15 by amylase, glycosidase and galactosidase;
- a kit according to the invention wherein the enzyme mixture furthermore comprises an enzyme chosen from the group formed by hydrolases that break ester bonds, or esterases, especially carboxyl ester
- 20 hydrolases, such as lipase, phospholipases and/or phosphonodiesterases, such as ribonuclease;
- a kit according to the invention, wherein the enzyme mixture is pancreatin;
- a kit according to the invention, wherein the
- 25 detergent exhibiting antibacterial, antiviral and antifungal properties is a solution containing surfactants;
- a kit according to the invention, wherein the detergent exhibiting antibacterial, antiviral and
- 30 antifungal properties is an alkaline solution containing surfactants and a quaternary ammonium;
- a kit according to the invention, wherein the detergent exhibiting antibacterial, antiviral and
- 35 antifungal properties is a neutral solution containing surfactants;
- a kit according to the invention, wherein the detergent exhibiting antibacterial, antiviral and
- antifungal properties is an acid solution containing surfactants;

- a kit according to the invention, which furthermore includes a solution of a disinfectant such as a sodium hypochlorite or potassium hypochlorite solution;

5 - a kit according to the invention, which furthermore includes a solution of an acid capable of dissolving deposits of mineral salts such as calcium carbonate;

10 - a kit according to the invention, wherein, in the acid solution for removing the deposits of mineral salts, the acid is chosen from the group formed by citric acid, peractetic acid, glycolic acid and oxyacetic acid.

15 In a variant of the invention, the solution comprising the enzyme mixture and the solution comprising the detergent form a single solution; the invention then relates to a composition intended for removing a biofilm, which method comprises an enzyme mixture
20 containing at least one enzyme chosen from the group of proteases and at least one enzyme chosen from the group of osidases, and a detergent exhibiting antibacterial, antiviral and antifungal properties.

25 In one particular embodiment, the invention relates to a composition in which the enzyme mixture is pancreatin.

Description of the preferred embodiments

30 The term "biofilm" is understood to mean a set of microorganisms that have grown on a support, especially bacteria, viruses, parasites and fungi. This biofilm grows and the microorganisms secrete a matrix of
35 exopolymers containing, inter alia, exopolysaccharides that will form a biological film called "slime" or "glycocalix" and is in the form of a gelatinous deposit on the surface of the walls.

The term "pancreatin" is understood to mean a pancreatic extract containing the digestive enzymes from the pancreas, especially proteolytic enzymes or proteases and enzymes such as lipase, amylase, ribonuclease and trypsin.

The term "detergent" is understood to mean any product whose composition has been specially designed to develop detergency phenomena, and which comprises, as essential components, surface agents, which are surfactants, and optionally additional components (various adjuvants, reinforcing agents, fillers and additives). Surfactants are chemical compounds which, when introduced into a liquid, lower its surface tension, with the effect of increasing the wetting properties.

The effectiveness of the method according to the invention was tested experimentally as described below.

The method was tested and implemented experimentally on five types of biofilm; biofilms 1, 2, 3 and 5 were obtained using an in vitro model that mimics the hemodialysis generator:

- biofilm 1: enriched with nutrients for accelerated growth (3 days) moderate thickness (about 10^5 CFU/cm²), very rich in slime;
- biofilm 2: not enriched with nutrients having grown over 1 month, equivalent to those actually encountered in dialysis generators (about 10^3 CFU/cm²), very rich in scale crystals;
- biofilm 3: enriched with nutrients for accelerated growth (5 days), thick (about 10^5 CFU/cm²), very rich in slime;
- biofilm 4: specimen of tubing conveying water for hemodialysis, taken from a center, covered with a biofilm of about 10^3 CFU/cm² but having grown over more than one year; and

biofilm 5: enriched with nutrients for accelerated growth, having grown in a "preventive" model over three weeks.

5 Production of the in vitro model

10 A 250 ml reactor vessel was filled with a nonsterile dialysate, prepared by diluting sterile apyrogenic hemodialysis concentrates (Clearflex® from Bieffe Medical) with nonsterile osmosed water containing *Pseudomonas putida*, *Pseudomonas fluorescens* and *Flavimonas orizibitans*, produced continuously in the laboratory.

15 The contaminating medium was circulated in a closed circuit in a loop of silicon tubing 1.5 meters in length and 5 mm inside diameter with a flow rate of 500 ml/min by a peristaltique pump. All the tubing and the reactor vessel were sterilized beforehand in an autoclave at 121°C for 30 minutes. Thus, the dialysate naturally contaminated by the bacteria in the water was the only source of microorganisms.

25 The entire system was maintained at a temperature of 37°C by a hotplate on which the reactor vessel was placed.

30 In the case of biofilms 1, 3 and 5, the bacterial growth and consequently the development of the biofilm were accelerated by adding an LB culture medium diluted 50-fold, i.e. an LB culture medium diluted to 5-fold and topped up with a flow equal to $1/10^{\text{th}}$ of that of the dialysate. The dialysate and culture medium top-up flow rates, regulated by a peristaltic pump, were 5 and 35 0.5 ml/minute, respectively.

In the case of biofilm 5, the model was modified in the following manner:

A nonsterile dialysate enriched with culture medium was made to flow for four hours through silicon tubing segments connected together by polypropylene couplers. Every four hours, the tubings were disconnected and integrated into disinfection systems (see below). After treatment, the tubings were reconnected and the contaminating medium resumed circulation for four hours. In parallel, control tubings that had never undergone disinfection were distributed in the circuit. Each day, two hemodialysis sessions, each interrupted by a disinfection session were thus able to be carried out. Overnight and at the weekends, the system was stopped after the last disinfection and the tubings were kept empty at room temperature. The system operated until a mature biofilm had developed on the control tubings.

Products and combinations tested

Seventeen products belonging to six different families were tested. The list of these products is given in table I. These products were evaluated singly or in combinations. Thus, a complete screening of sixty combinations was carried out on biofilm 1; nine combinations were then evaluated on biofilm 2; finally, the best combination selected was tested on biofilms 3, 4 and 5.

Table I: List of Products Tested

Family	Product	Supplier
Surfactants/detergents	Sodium dodecyl sulfate	Sigma
	Triton	Sigma
	RBS	Chemical products
	Tween	Sigma
Enzymes	Trypsin	Sigma
	Pancreatin	Sigma
	Fungal protease	Sigma
	Thermolysin	Sigma
Acids	Perchloric acid	Merck
	Citric acid	Merck
	Trichloroacetic acid	Merck
Cell dissociation products	Versene	Sigma
	Cell dissociation	Sigma
Alkalis	NaOH	Prolabo
	KOH	Prolabo
Miscellaneous	Bleach	-
	pH10 buffer (bicarbonate)	Prepared in the lab

Sampling

- 5 The tubings covered with biofilms 1, 2 and 3 were cut into segments 5 cm in length. For screening on biofilms 1 and 2, each segment for undergoing one of the various treatments to be investigated was selected by drawing
- 10 lots. Control specimens taken at random from the silicon loop were kept untreated.

Treatment of the specimens

- 15 The tubing segments to be treated were attached to the descending branch of a set-up consisting of two tubings - one ascending and the other descending, a peristaltic pump and a water bath (for treatments carried out a temperature above 20°C). The product to
- 20 be tested in "recirculation" mode was dissolved in a 100 ml flask and driven by the peristaltic pump at a rate of 500 ml/min in a closed circuit through the

tubings for a duration corresponding to the contact times given in table II. The product to be tested in "static" mode was dissolved in a 100 ml flask and driven by the peristaltic pump until the tubings had been filled; the pump was then stopped and the product left in stasis for the desired contact time. After each treatment by a given product, the tubing specimens were rinsed for five minutes with osmosed water.

10 Methods of evaluating the effectiveness of the treatments:

Three fundamental parameters were used to evaluate the effectiveness of the treatments:

- 15 The reduction in area covered;
 The reduction in number of culturable bacteria;
and
 The reduction in the level of endotoxins.

- 20 The screening on biofilms 1 and 2 only took the first parameter into account. The best combination adopted was then evaluated in depth as regards its effectiveness on bacterial mortality and endotoxin elimination.

25

Method of quantifying the area covered:

The control and treated biofilms were stained:

- 30 - either with a 0.25% crystal violet solution;
 - or with a BacLight® fluorochrome solution (Syto 9 and propidium iodide). The viable bacteria show up green while the dead bacteria show up yellow or red.

- 35 The silicon tubing specimens covered with stained biofilms were attached to glass slides and observed under an optical microscope, the microscope being connected to a camera and to "Scion Images" image analysis software. Thus, several (6 to 10) photographs of the same specimen were taken, the stained area was

evaluated quantitatively by the image analysis software and a mean covered area value per specimen was calculated. This mean value was compared with the covered area of the untreated control specimens: a
5 percentage reduction in covered area was then calculated.

Moreover, for more accurate examination, the specimens treated by the most effective combination were observed
10 under a laser confocal microscope.

Method of quantifying the culturable bacteria:

The biofilm covering the tubing specimens was detached
15 from the substrate using a mechanical scraper, ensuring complete, uniform and reproducible detachment. This scraper consisted of a power screwdriver at the end of which a flame-sterilizable stainless steel spatula was fixed. By rotating the spatula in the lumen of the
20 tubing, the biomass was driven to the bottom of a sterile tube. This action was facilitated by a stream of sterile water. Any bacterial aggregates were then separated via the needle of a syringe. The number of culturable bacteria was determined by counting the CFUs
25 after plating the resulting bacterial suspension on R₂A agar and incubating at room temperature for 7 days.

More precisely, the specimens proving to be uncontaminated after plating were completely filtered
30 and the filtration membrane was incubated on R₂An agar at room temperature for 7 days.

Endotoxin assay method:

35 The bacterial endotoxins were quantified in the bacterial suspension resulting from the detachment (see above) by a standardized reference test, namely the kinetic chromogenic LAL test (Charles River Endosafe).

Results of the study on biofilms 1 and 2:

The results of the screening on biofilms 1 and 2 are shown in tables II and III.

5

Table II: Screening on biofilm 1

Treat. No.	Product	Conc.	°C	Time	Mode	Result
1	SDS	5%	RT	40 min	Stat	**
2	Triton	5%	RT	40 min	Stat	*
3	RBS	5%	RT	40 min	Stat	**
4	Tween	5%	RT	40 min	Stat	.
5	Versene	Pure	RT	40 min	Stat	.
6	Trypsin	EDTA IX	RT	40 min	Stat	.
7	Trypsin	0.25%	37	40 min	Stat	**
8	SDS	5%	RT	40 min	Recirc	*
9	Tween	5%	RT	40 min	Recirc	.
10	RBS	5%	RT	40 min	Recirc	***
11	Perchlo ac	0.05%	RT	40 min	Stat	.
12	NaOH	0.01N	RT	40 min	Stat	*
13	KOH	0.02N	RT	40 min	Stat	**
14	TCA	0.25%	RT	40 min	Stat	.
15	KOH	0.02N	RT	40 min	Recirc	**
16	Triton	5%	RT	40 min	Recirc	.
17	RBS	5%	RT	1 h	Recirc	***
18	RBS	5%	RT	24 h	Recirc	****
19	RBS	2%	RT	1 h	Recirc	***
20	RBS	2%	RT	24 h	Recirc	****
21	KOH	0.02N	RT	1 h	Recirc	***
22	KOH	0.02N	RT	24 h	Recirc	****
23	Pancreatin	1%	37	1 h	Recirc	****

24	Pancreatin	1%	37	1 h	Stat	**
25	Pancreatin	0.10%	37	2 h	Recirc	.
26	Citric acid	3%	RT	40 min	Recirc	.
27	KOH	0.002N	RT	40 min	Recirc	**
28	Cell dissoc	Pure	37	40 min	Stat	**
29	Trypsin	0.25%	37	40 min	Recirc	.
30	Protease	0.25%	37	40 min	Recirc	**
31	RBS	2%	RT	40 min	Recirc	****
32	RBS + Cl	2%+0.2%	RT	40 min	Recirc	****
33	Thermolys	2mg/50ml	37	40 min	Recirc	**
34	Pancreatin	0.50%	37	40 min	Recirc	***
35	KOH	0.001N	RT	40 min	Recirc	**
36	RBS	2% pH7	RT	40 min	Recirc	***
37	Trypsin	1%	37	40 min	Recirc	*
38	Protease	1%	37	40 min	Recirc	**
39	RBS	1% pH10	RT	40 min	Recirc	****
40	RBS	0.5% pH10	RT	40 min	Recirc	***(*)
41	RBS	1% pH10	RT	5 min	Recirc	***
42	RBS	0.1% pH10	RT	40 min	Recirc	**(*)
43	RBS + Pancreatin	1% pH10 1%	40	5 min	Recirc	*
44	RBS + Pancreatin	1% pH10 0.50%	40	5 min	Recirc	**
45	RBS + Pancreatin	1% pH10 0.50%	40	40	Recirc	***
46	pH10 Buffer	Pure	37	40	Recirc	*
47	Pancreatin	0.5pH10	37	40	Recirc	**
48	Pancreatin and RBS	0.5% pH7 1% pH10	37 RT	5 min 5 min	Recirc Recirc	*****
49	Pancreatin	0.25% pH7	37	5 min	Recirc	***
50	Pancreatin + Thermolys	0.25% pH7 2mg/50ml	37	5 min	Recirc	**

51	Pancreatin + Thermolys + Protease	0.25% pH7 2mg/50ml 0.25%	37	5 min	Recirc	**
52	Pancreatin + Thermolys + Protease + Trypsin	0.25% pH7 2mg/50ml 0.25% 0.25%	37	5 min	Recirc	**
53	Pancreatin and RBS	0.25 1%	37 RT	5 min 5 min	Recirc	*****
54	Pancreatin and RBS	0.25 0.50%	37 RT	5 min 5 min	Recirc	*****
55	Pancreatin and RBS	0.25% 0.50%	37 RT	5 min 30 min	Recirc	*****
56	Pancreatin and RBS	0.10% 0.50%	37 RT	5 min 5 min	Recirc	*****
57	Pancreatin and RBS	0.10% 0.10%	37 RT	5 min 5 min	Recirc	*****
58	Citric acid and RBS	3% 1%	37 RT	5 min 30 min	Recirc	*****
59	Citric acid and pancreatin and RBS	3% 0.50% 1%	37 37 RT	5 min 5 min 30 min	Recirc	*****
60	Pancreatin and RBS	0.50% 1%	37 RT	5 min 30 min	Recirc	*****

Key	Elimination	% covered area reduction
.	No removal	0
*	Poor removal	0-25%
**	Moderate removal	25-50%
***	Good removal	50-75%
	Excellent removal	75-99%
*****	Complete removal	100%

+ = mixture;

and = sequential application

Table III: Screening on biofilm 2

Treat. No.	Product	Conc.	°C	Time	Mode	Result
1	Pancreatin	0.5% pH 7	37	5 min	Recirc	***
	and RBS	1% pH10	RT	5 min	Recirc	
2	Citric acid and RBS	3%	37	5 min	Recirc	*****
		1%	RT	30 min		
3	Pancreatin and citric acid and RBS	0.50%	37	5 min	Recirc	*****
		3%	37	5 min		
		1%	RT	30 min		
4	Pancreatin and RBS	0.50%	37	5 min	Recirc	***
		1%	RT	30 min		
5	Citric acid and pancreatin and RBS	3%	37	5 min	Recirc	*****
		0.50%	37	5 min		
		1%	RT	30 min		
6	Pancreatin and RBS	0.50%	37	5 min	Recirc	***
		1%	RT	30 min		
7	Citric acid	3%	37	5 min	Recirc	****
8	Citric acid and RBS	3%	37	5 min	Recirc	*****
		0.50%	RT	30 min		
9	Citric acid and RBS	3%	37	5 min	Recirc	****
		0.10%	RT	30 min		

5 The combination adopted following these screening tests was that which gave the best removal of the two biofilms, namely the following "combination K":

10 Product A = pancreatin®, a laboratory reactant sold by Sigma: Pig pancreas extract, namely an enzyme mixture containing, inter alia, lipase, protease, amylase, trypsin, ribonuclease, etc.

Product B = citric acid

15 In the particular case of disinfecting dialysis generators, this product acts as a decalcifier and removes the scale crystals that trap the bacteria and promote adhesion of the biofilm to the substrate.

Product C = RBS®, a foaming alkaline detergent solution, sold by Chemical Products, exhibiting bactericidal, virucidal and fungicidal properties, containing surfactants and a quaternary ammonium (disinfectant).

Qualitative and quantitative data:

- 10 Photographs of biofilms 1 and 2 before and after the action of combination K have allowed the action of the combination to be visually quantified.

15 Table IV gives the values of the parameters measured before and after action of combination K.

Tables IV: Quantitative data for evaluating the effectiveness of combination K on biofilms 1 and 2

20 Table IVa): Biofilm 1

Parameter	Before treatment	After treatment	% reduction
Area covered (sq.in)	20	<0.001	>99.99
Culturable bacteria (CFU/cm ²)	10 ⁵	<1	>99.999
Endotoxins (EU/cm ²)	10039	<0.005	>99.99

Table IVb): Biofilm 2

Parameter	Before treatment	After treatment	% r duction
Area covered (sq.in)	13.4	<0.001	>99.99
Culturable bacteria (CFU/cm ²)	3x10 ³	<1	>99.999
Endotoxins (EU/cm ²)	40	<0.005	>99.99

Determination of the MIC of RBS:

5

Since the anitbacterial capacity of this combination is provided by RBS, its minimal inhibitory concentration was determined on the microorganisms constituting the biofilms under examination.

10

A mixture of contaminated fresh dialysate (prepared with nonsterile osmosed water containing the microorganisms described above) and LB medium in the proportions of 50/50 v/v was prepared. RBS solutions with concentrations of 100%, 50%, 10%, 5%, 1%, 0.5% and 0.1% were produced by cascade dilutions, and then 300 µl of each of these solutions were added to 3 ml of the contaminated mixture. After incubation for 12 hours at room temperature, the CFUs were counted on R₂An agar for each of the RBS concentrations tested.

20

The MIC is defined as being the lowest concentration that inhibits the growth of microorganisms. The results are given in table V.

25

Table V: Determination of the MIC of the RBS

Conc (%)	0	0.01	0.05	0.5	1	2	3	4	5	10
UPC/ml	2x10 ⁸	4.8x10 ⁷	2.1x10 ⁷	9x10 ⁶	6.1x10 ⁵	1.4x10 ⁵	1200	200	0	0

For safety, the choice was made to use the RBS solution in dilute form so as to obtain an MIC of 1.5.

Initial protocol adopted:

5

0.5%/ pH 7.3 pancreatin: preparation for 100 ml: 500 mg of pancreatin powder + 1 g of powdered PBS (phosphate buffer saline) buffer (Sigma), diluted in 100 ml of Hemodialysis Water (HDW).

10

Closed-circuit flow of the solution in the tubings with a flow rate of 500 ml/min for 5 minutes at 37°C.

15 Rinsing: 5 minutes with HDW at 500 ml/min in open circuit.

3%/pH 2.2 citric acid: preparation for 100 ml: 3 g of powdered citric acid (Merck) in 100 ml of HDW.

20

Closed-circuit flow of the solution in the tubings with a flow rate of 500 ml/min for 5 minutes at 20°C.

Rinsing: 5 minutes with HDW at 500 ml/min in open circuit.

25

7%/pH 10 RBS®: preparation for 100 ml: 7 ml of concentrated solution + 566 mg of powdered sodium carbonate + 388 mg of powdered sodium bicarbonate, diluted with HDW to make up to 100 ml.

30 Closed-circuit flow of the solution in the tubings with a flow rate of 500 ml/min for 30 minutes at 20°C.

Rinsing: 5 minutes with HDW at 500 ml/min in open circuit.

Results of the study of biofilms 3 and 4:

35

Combination K in its initial form described above left a few cells adhering to biofilms 3 and 4, particularly thick or old biofilms. To remove such biofilms, an "enriched formula" of combination K was developed.

"Enriched" formula:

1 1%/pH 7.3 pancreatin: preparation for 100 ml: 1 g of
5 pancreatin powder + 1 g of powdered PBS (phosphate
buffer saline) buffer (Sigma), diluted in 100 ml of
hemodialysis water (HDW):

10 Closed-circuit flow of the solution in the tubings
with a flow rate of 500 ml/min for 5 minutes at 37°C.

Rinsing: 5 minutes with HDW at 500 ml/min in open
circuit.

15 5%/pH 2.2 citric acid: preparation for 100 ml: 5 g
of powdered citric acid (Merck) in 100 ml of HDW.

Closed-circuit flow of the solution in the tubings with
a flow rate of 500 ml/min for 30 minutes at 20°C.

20 Rinsing: 5 minutes with HDW at 500 ml/min in open
circuit.

25 15%/pH 10 RBS®: preparation for 100 ml: 15 ml of
concentrated solution plus 566 mg of powdered sodium
carbonate + 388 mg of powdered sodium bicarbonate
+ 6 ml of 5.2% concentration bleach (final sodium
hypochlorite concentration: 0.3%), diluted in HDW to
make up to 100 ml.

30 Closed-circuit flow of the solution in the tubings
with a flow rate of 500 ml/min overnight at 20°C.

Rinsing: 5 minutes with HDW at 500 ml/min in open
circuit.

Qualitative and quantitative data:

Photographs of biofilm 4 before and after the action of enriched combination K allowed the effectiveness of the method according to the invention to be visually quantified.

Table VI gives the values of the parameters measured before and after the action of enriched combination K on biofilm 4.

Table VI: Quantitative data for evaluating the effectiveness of enriched combination K on biofilm 4

Parameter	Before treatment	After treatment	% reduction
Area covered (sq.in)	25	<0.001	>99.99
Culturable bacteria (CFU/cm ²)	1600	<1	>99.999
Endotoxins (BU/cm ²)	115	<0.005	>99.999

Results of the study on biofilm 5:

A very thick (more than 3×10^8 CFU/cm²) biofilm very rich in slime, very rich in bacterial endotoxins and completely covering the surface of the specimen (30 square inches) grew on the surface of the untreated control specimens, whereas only a few adherent dead cells were deposited on the surface of the treated specimens, treated every 4 hours with unenriched combination K (initial formula).

The quantitative data are given in table VII.

Table VII: Effectiveness of combination K on biofilm 5

Parameter	Without treatment	With treatment	% inhibition
Area covered (sq.in)	30	1.3	96
Culturable bacteria (CFU/cm ²)	3.9x10 ⁹	<1	>99.999
Endotoxins (EU/cm ²)	65282	0.4	>99.999

5 Photographs of the control biofilm and of the treated specimens allowed the effectiveness of the method according to the invention to be verified.

10 The method and the kit according to the invention can be used in the circuits of hemodialysis equipment to combat legionellosis for example in hot-water circuits and air-conditioning systems and cooling towers, in the agrifoodstuffs industry, in climate-controlled rooms or in confined-atmosphere rooms, for cleaning dentistry equipment, and for reusable and non-autoclavable
15 medical instruments.

In fluid flow equipment, the method according to the invention will be employed by introducing the solution or solutions simultaneously or sequentially into the
20 circuits from which the biofilm must be removed, by making said solution(s) circulate for a period sufficient to allow the biofilm to be removed, followed by purging and rinsing if necessary.

25 For treating surfaces, worksurfaces and prostheses, the method according to the invention will be carried out by application of or by immersion in the solution or solutions according to the invention, sequentially or simultaneously, followed if necessary by rinsing.

**This Page is Inserted by IFW Indexing and Scanning
Operations and is not part of the Official Record**

BEST AVAILABLE IMAGES

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images include but are not limited to the items checked:

- ☐ BLACK BORDERS
- ☐ IMAGE CUT OFF AT TOP, BOTTOM OR SIDES
- ☐ FADED TEXT OR DRAWING
- ☐ BLURRED OR ILLEGIBLE TEXT OR DRAWING
- ☒ SKEWED/SLANTED IMAGES
- ☐ COLOR OR BLACK AND WHITE PHOTOGRAPHS
- ☐ GRAY SCALE DOCUMENTS
- ☒ LINES OR MARKS ON ORIGINAL DOCUMENT
- ☐ REFERENCE(S) OR EXHIBIT(S) SUBMITTED ARE POOR QUALITY
- ☐ OTHER: _____

IMAGES ARE BEST AVAILABLE COPY.

As rescanning these documents will not correct the image problems checked, please do not report these problems to the IFW Image Problem Mailbox.